era**visc** X

VISCOSITY & DENSITY TESTING

eravisc X Xtra Easy Xtra Rugged Xtra Accurate

eravisc X is the world's smallest and most rugged kinematic viscometer with lab-grade precision meeting and ASTM D7042. ASTM D445 Its revolutionary viscosity cell is combined with eralytics unique x-oscillating metal U-tube density cell. This allows an independent highprecision density measurement in full compliance with ASTM D4052 and ISO 12185.



eravisc X's revolutionary capillary measuring cell is a real game changer for high-precision kinematic viscosity testing at any temperature between 15 °C and 100 °C. Latest Peltier Technology allows for the automated, high-speed viscosity index measurement. eralytics unique FillingProof™ solution detects even smallest bubbles, ensuring maximum density accuracy.

eravisc X

Sample

Rinse

Method Viscosity

Operator eralytics

Lube 1

T Dens [°C 15.00

eralytics 110

40.00 °C

15.00

eravisc X is the analyzer with the smallest footprint on the market, still offering maximum ease of use and operator convenience. The compact dimensions and lightweight instrument design offer true portability on-site and require minimum bench space in the lab.

The fully automated measurement of the kinematic viscosity at 40 °C or at 100 °C takes only 1 minute. Results are displayed on the large color touchscreen and can be printed, exported via USB or directly transferred to a PC or LIMS via Ethernet connection.

With its unique sample introduction from the top eravisc X is equally easy to operate for right-handed and left-handed persons.

eravisc X is built in a corrosion-resistant aluminum housing. Its robustness and reliability make it the ideal solution for daily viscosity measurements in the laboratory as well as directly in the field to immediately prove specification or detect contamination, and oil degradation.

eravisc X is the perfect solution for highly accurate testing even under extreme conditions.

> dependent measurements of viscosity and density enabling the simultaneous determination of both

eralytics unique, dynamic capillary viscosity measurement method determines the flow time of the sample as a function of a rapidly measured pressure. This patented principle combines the advantages of traditional capillary methods with the speed and ease of use of eralytics instruments.

eravisc X allows for the invalues at different temperatures.

Technical Specifications of era**visc X**

)
)
)
193,417 B2)
intenance
B dongle

^{*} Typical values, depending on sample type and temperature and without consideration of the particular uncertainties of reference materials (CRM).

Due to continuing product development, specifications are subject to change.

All eralytics products are manufactured under ISO 9001 regulations and are CE, ROHS and UL/CSA compliant. www.eralytics.com/eravisc-x



eralytics instruments are available worldwide.

An international network of over 50 authorized and well-trained distributors is ready to answer your inquiries and to offer local support and service.

www.eralytics.com/distribution

eralytics⁰

Autokaderstrasse 29, Building 4A 1210 Vienna, Austria Phone: +43 1 890 50 33 0 office@eralytics.com www.eralytics.com